

YDFIX Distal Femur



Prof.Dr. Yunus DOGRAMACI

M.Kemal University

Prof. Dr. Aydın KALACI

M.Kemal University

Ass.Prof.Dr. Ali BİLGE


Canakkale University

1.	Introduction	4	4.	Implant list	16-17
1.1	The implant	4	5.	Instrument list	18
1.2	The instruments	4			
1.3	Indications	4			
2.	Implant range	5			
3.	Surgical description	6-9			
3.1	Patient positioning	6			
3.2	Plate selection	6			
3.3	Assembly of plate and targeting arm	6			
3.4	Incision	6			
3.5	Modellation of the plate	7			
3.6	Insertion of the plate	7			
3.7	Positioning with the oval hole	7			
3.8	Closing the frame	9			
3.9	Temporary fixation of the plate	9			
3.10	Fixation of the head	9			
3.11	Planned direction locking	10			
3.11.1	Assembly of the head block	10			
3.11.2	Planned direction drilling – head	10			
3.11.3	Length gauging – head	10			
3.11.4	Screw insertion – head	11			
3.11.5	Drilling – tail	11			
3.11.6	Length gauging – tail	12			
3.11.7	Screw insertion – tail	12			
3.12	Selected direction locking	13			
3.12.1	Selected direction drilling – head	13			
3.12.2	Length gauging – head	13			
3.12.3	Screw insertion – head	14			
3.12.4	Drilling – tail	14			
3.12.5	Length gauging – tail	15			
3.12.6	Screw insertion – tail	15			

Ydfix plates, serving to heal peri- and intraarticular fractures, have a new family member, **YDF (Ydfix Distal Femur)** which goes to the distal part of the femur bone offering the same usual high end features that we got used to from Ydfix plates. The comfort and security of polyaxial locking is further enhanced by a minimally invasive radiolucent targeting arm for the tail. To reach stable-enough locking the system uses dia 5,1 mm polyaxial and cortical screws. They represent the outstanding features of Ydfix screws together with excellent mechanical properties.

1.1 | The implant

- Polyaxial angle stabilized system in step - free ± 15 deg angulation of insertion
- Optimal, pre-determined screw directions in the holes
- Maximum 3 times of correction possibility when misidentifying the correct screw direction
- Thinned head, the implant does not interfere with the soft tissues
- Rounded edges to protect nearby soft tissues
- Oval hole for plate positioning
- Ability to perform minimally invasive surgery
- Self tapping but blunt ended screws to avoid tissue irritation



- Anodized Titanium raw material
- Torx recess screws

1.2 | The instruments

- Capable of drilling in preset and ± 15 deg directions step - free
- Easy to assemble targeting head for choosing either selected or planned direction fixation in the pre-determined directions.
- Radio translucent targeting arm for the tail holes
- Instruments and implants in one tray
- Optimized instruments
- Color coded torque limiting screwdriver

1.3 | Indications

Distal femoral fractures
Periprosthetic fractures

2.1 | Ydfix Distal Femur plate

Holes on tail	Side
6H	right/left
8H	right/left
10H	right/left
12H	right/left
14H	right/left
16H	right/left



Raw material
anodized Titanium

Color
grey

2.3 | Ydfix screw Ø5,1 mm



Length (mm)
24 - 80

Raw material
anodized Titanium

Color
purple

2.4 | Ydfix screw Ø6,5 mm



Length (mm)
25 - 80

Raw material
anodized titanium

Color
grey

2.2 | Ydfix screw Ø5,1 mm



Length (mm)
12 - 20

Raw material
anodized Titanium

Color
purple

2.4 | Cortical screw - TX Ø5,1 mm



Length (mm)
24 - 55

Raw material
anodized Titanium

Color
grey

3 | Surgical description

3.1 | Patient positioning

In supine position on straight, translucent table with the possibility of bi-lateral x-ray control.

3.2 | Plate selection

During preoperative planning select the most appropriate plate for the fracture and the anatomy of the patient. Consider that a little longer plate is more acceptable than a too short one.

3.3 | Assembly of plate and targeting arm

The system enables selected- and planned direction locking on the head of the plate. See point 3.11.1 for the planned direction head targeting.

Place the targeting arm on the plate while paying attention to the small peg entering into its hole on the plate. Mount the carbon arm and fix the system to the plate as per the image.

Check accuracy with soft tissue protector, drill sleeve and the dia 4 mm drillbit.

In case we do not use the longest plate put the indicator plug into the arm to represent the position of the last hole.



3.4 | Approach

Concerning the approach and the incision, the number of screws applied and the optimal steps of the surgery this present description does not make a stance. The above shall be acquired from surgical textbooks, videos and workshops.

3.5 | Modellation of the plate

For optimal fitting we can modellate the plate if needed. Use table bending device for that purpose.

Attention!

After bending the plate the targeting devices on the head and the tail cannot be used!

3.6 | Insertion of the plate

The plate can be introduced minimally invasively. Guide the plate mounted on the targeting arm on the bone surface towards the proximal. Optimal position shall be checked with image intensifier.

The plate can be fixed with Kirschner wires at this step but the fine tuning of the position with the oval hole is only possible after the wires are removed.



3.7 | Positioning with the oval hole

Lock the plate in the oval hole with a grey cortical screw.


Pre drilling takes place through the straight side of the Ø4 mm double drill sleeve. After setting the optimal angle perform drilling with the Ø4 mm drillbit.

3 | Surgical description

Measure the necessary screw length.

Remove the Ø4 mm double drill sleeve. Hook the gauge to the other side of the hole while moving the reader on the bone surface. Read length at the red mark.

The plate is fixed to the bone with a Ø5,1 mm cortical screw in the oval hole. The screw is not fully tightened until the fine tuning of the position is done. Afterwards it is to be locked firmly.




3.8 | Closing the frame

Perform incision at the most proximal hole of the plate used. Push the Ø8 mm soft tissue protector into the hole and turn the drill sleeve into the hole of the plate. This way the plate and the carbon arm are aligned. The plate can be fixed - either at the most proximal hole or through the small hole on the tip of the plate - with Kirschner wires to the bone.

Attention!


This step is absolutely necessary for the perfect targeting.

In case sleeves cannot be turned by the hand, use the T25 screwdriver!



3.9 | Temporary fixation of the plate

Fix the plate to the bone with Kirschner wires through the corresponding holes on the head or at the tip of the plate. In the latter the wires can be led from the sides of the carbon arm.



3.10 | Fixation of the head

There are 6 locking options on the head. In each hole planned- and selected direction locking is also possible.

3.11 | Planned direction locking

In this case the screws can be inserted in the anatomically optimal direction. The screw-plate connection will be angle stabilized ensuring loosening-free locking.


3 | Surgical description

3.11.1 | Assembly of the head block

Fix the removable head as per the image.


The mounted targeting arm for the tail ensures the targeting head against rotation. The targeting head and the plate are connected with a short sleeve in the center hole.

You can use the double drill sleeve's straight side for monoaxial locking as well. This technique is to be detailed further on.



3.11.2 | Planned direction drilling - head

Drill through the sleeve for the Ø5,1 mm screw while using image intensifier control. The spiral drill to be used is Ø4 mm in diameter.




3.11.3 | Depth gauging – head

Depth gauging can take place in two ways.

As per the first, previously already detailed, use a hooked gauge. Hook it into the other wall of the hole, push the moving part on the bone surface and read length.




The other method uses the purple drill stop over the Ø4 mm drillbit just above the spiral part. Perform drilling through the sleeve. The necessary length can be read on the scale of the drillbit at the drill stop's side facing the sleeve.



3.11.4 | Screw insertion – head

Drive the selected Ø5,1 mm screw with the T25/3,5 Nm torque screwdriver. Using torque screwdriver reduces chances of cold welding during the healing period.

Repeat the above steps in case of all necessary screws. This way they are inserted in the anatomically optimal direction.



3.11.5 | Drilling – tail

YDF system uses radio translucent arm for minimally invasive targeting of the tail.


Attention!

The targeting arm can be used only with intact, non-modellated plate. Otherwise drilling accuracy greatly reduces and plate and/or targeting arm damages might occur.

3 | Surgical description


On the tail the traditional planned direction locking is suggested.

Push and drive the soft tissue protector and the drill sleeve into the hole's thread you wish to use. Perform drilling with the Ø4 mm drillbit. Remove the Ø4 mm drill sleeve.




3.11.6 | Depth gauging – tail

To determine screw length either use depth gauge or the drill stop.



3.11.7 | Screw insertion – tail

Turn in the selected screw with the T25/3,5 Nm screwdriver.




3.12 | Selected direction locking

In case of selected direction locking we have the freedom of ± 15 degrees from the anatomically optimal direction. Use the system without the removable head to have the possibility of the six selected direction locking on the head.

3.12.1 | Selected direction drilling – head

Place the conical end of the Ø4 mm double drill sleeve into the hole. It fits exactly into the hole and its symmetry axis is in the anatomically optimal direction (the same direction in which monoaxial locking takes place). Thus it is ensured that the ± 15 degree direction is always from the anatomically optimal direction.

Perform drilling in the desired direction with the Ø4 mm drillbit.




3.12.2 | Length gauging – head

Use the depth gauge without the sleeve.

Attention!


The drillstop method cannot be used when drilling through the conical end of the double drill sleeve.



3 | Surgical description

3.12.3 | Screw insertion – head

Drive in the selected screw with the T25/3,5 Nm torque limiting screwdriver.



3.12.4 | Drilling – tail


YDF system uses radio translucent arm for minimally invasive targeting of the tail.

Attention!

The targeting arm can be used only with intact, non-modellated plate. Otherwise drilling accuracy greatly reduces and plate and/or targeting arm damages might occur.


On the tail the traditional monoaxial locking is suggested.

Push and drive the soft tissue protector and the drill sleeve into the hole's thread you wish to use. Perform drilling with the Ø4 mm drillbit. Remove the Ø4 mm drill sleeve.




3.12.5 | Length gauging – tail

To determine screw length either use depth gauge or the drill stop.



3.12.6 | Screw insertion – tail

Turn in the selected screw with the T25/3,5 Nm screwdriver.




Compression locking

For implanting compression screws drill through the compression side of the double drill sleeve V with the 4 mm drillbit. The arrow on the sleeve shall point to the direction of the fracture. Take general guidelines of compression locking into account. Length gauging is followed by screw insertion with the T25 screwdriver.

For the optimal compression use TX cortical screws.

Gauging tip

When using the drill stop gauging the 4x250 mm bit is calibrated to the short sleeves while the 4x260 mm bit is for the long sleeves. Due to scale positions no false reading is possible.



4.1 | Ydfix Distal Femur plate



Anodized Titanium

Cat no	Size
3011402006	6H/left
3011402008	8H/left
3011402010	10H/left
3011402012	12H/left
3011402014	14H/left
3011402016	16H/left
3011401006	6H/right
3011401008	8H/right
3011401010	10H/right
3011401012	12H/right
3011401014	14H/right
3011401016	16H/right

4.2 | Ydfix screw Ø5,1 mm



Anodized Titanium

Cat. no	Size (mm)
1017451024	24
1017451026	26
1017451028	28
1017451030	30
1017451032	32
1017451034	34
1017451036	36
1017451038	38
1017451040	40
1017451042	42
1017451044	44
1017451046	46
1017451048	48
1017451050	50
1017451055	55
1017451060	60
1017451065	65
1017451070	70
1017451075	75
1017451080	80

4.3 | Ydfix screw - blunt Ø5,1 mm

Anodized Titanium

Cat. no	Size (mm)
1017451012	12
1017451020	20



4.5 | Cortical screw - TX Ø5,1

Anodized Titanium

Cat. no	Size (mm)
1032451024	24
1032451026	26
1032451028	28
1032451030	30
1032451032	32
1032451034	34
1032451036	36
1032451038	38
1032451040	40
1032451042	42
1032451044	44
1032451046	46




4.4 | Ydfix screw Ø6,5 mm

Anodized Titanium



Cat no	Size (mm)
1057465065	65
1057465070	70
1057465075	75
1057465080	80



5.1 | Instruments

Screwdriver (T25)	1 pc	 5210720025
Torque screwdriver (T25/3,5 Nm)	1 pc	 5210510046
Spiral drill (4x250 mm; 4x260 mm)	1-1 pc	 5280251903; 5280122915
Double drill sleeve - PAS (4 mm)	1 pc	 5280122910
Double drill sleeve - V (Large)	1 pc	 5275214902
Kirschner wire (2x230 mm)	5 pcs	 5937520230
Screw forceps	1 pc	 5939999002
Drill stop (4 mm)	2 pcs	 5210510240
Depth gauge (10-90 mm)	1 pc	 5280122912
Target device - YDF	1 pc	 5280251902
Filled up (YDF)		5233800024

Optional

Power torque screwdriver (T25/3,5 Nm) (for purple screw)	 5210510048
Spiral drill with quick connecting end (3,2x195 mm)	 5939532195
Screwdriver (3,5 mm)	 5210700035
Double drill sleeve (6,5/3,2 mm)	 5233910002



YDFIX

Isiso San Sit 13. Yol V2 No:20
Esenyurt-İstanbul / TURKEY
Phone: +90 212 438 63 67
Email: info@imd.com.tr



www.imd.com.tr