


# PROXIMAL FEMUR



# Why Ydfix Proximal Femur plate?

- ✓ Can be combined with cerclage system.
- ✓ Fits to the well-known Ydfix system both in terms of screws and instruments.
- ✓ Gives a  $\pm 15$  degree freedom to the surgeon, in total 30 degrees of cone angle for polyaxial locking.
- ✓ Plate shape is anatomic, the arched head fits to the trochanter tip.
- ✓ Protects soft tissues: self tapping but blunt screws, rounded contours.
- ✓ There are combined holes, the holes on the tail are also polyaxial.
- ✓ Supports minimally invasive technique.



İsıso San Sit 13. Yol V2 No:20  
Esenyurt-İstanbul / TURKEY  
Phone: +90 212 438 63 67  
Email: [info@imd.com.tr](mailto:info@imd.com.tr)